WOUNDCLOT™

Advanced Bleeding Control™

OPLOSBAAR HEMOSTATISCH GAAS

Steriele inhoud

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SOLUBLE HEMOSTATIC GAUZE

Sterile contents

**Description:** WoundClot ™ is a soluble gauze, based on cellulose structure. It is fully biocompatible and provided sterile. In contact with liquids it forms an adhesive gel which expands and adheres to the wound for bleeding control.

WOUNDCLOT™

Hemostatic Gauze

Advanced Bleeding Control™

INSTRUCTIONS FOR USE (CE)

FOR TEMPORARY CONTROL OF BLEEDING

Sterile contents

**Product size:** 2.5x2.5cm (1x1"), 5x5cm (2x2"), 5x10cm (2x4"), 10x10cm (4x4"),Twin pack- 2.5x2.5cm (1x1") & 5x5cm (2x2"). **Indications:** WoundClot™ gauze is intended for external bleeding control of superficial wounds, minor cuts, and abrasions. **Application:** Place the gauze onto the wound, apply slight pressure and let the blood flow through the gauze until it forms into gel and blood stops flowing. If needed, you may slightly wet the gauze to create a complete gel membrane. **Removal:** Once hemostasis is achieved, remove by saturating the gauze with water or saline solution, then gently wipe any excess gauze material away from the wound. **Cautions:** It is advised to clean the wound with water before using WoundClot™. In case of itching, burning, rash or signs of infection such as fever, pain, redness or swelling, seek professional medical care.

Do not re-sterilize

WoundClot™ gauze is a single use device. Do not reuse.

Consult instruction for Use

Sterilization using Irradiation

Single sterile barrier system with protective packaging inside

Do not use if package is damaged. Sterility of contents guaranteed unless individual package is damaged or opened.

Do not use after expiration date.

Store in dry place

Batch code

Catalog number

Medical Device

Unique Device Identifier

Manufacturer

CORE SCIENTIFIC CREATIONS

0483

Manufacturer: Core Scientific Creations Ltd.  
Office: 29 Nahshon Street, Petah Tikva 4926918, Israel  
Manufacturing Site: 6 Bahat Street, Petah Tikva, 4927149, Israel.  
T: +972-3-7359555 | F: +972-3-7359544  
info@wound-clot.com | www.woundclot.org

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MedNet EC-REP III GmbH
MedNet Swiss GmbH