WOUNDCLOT™

Advanced Bleeding Control™

OPLOSBAAR HEMOSTATISCH GAAS

Steriele inhoud

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SOLUBLE HEMOSTATIC GAUZE

Sterile contents

Description: WoundClot™ is a water soluble gauze, based on cellulose structure. It is fully biocompatible and provided sterile. In contact with liquids it forms an adhesive gel which expands and adheres to the wound for bleeding control.

English

WOUNDCLOT™

Surgical

INSTRUCTIONS FOR USE

SURGICAL DEVICE (CE)

Advanced Bleeding Control™

SURGICAL HEMOSTATIC DRESSING

For surgical procedures

Sterile contents

Product size: 2.5X2.5cm (1X1”), 5X5cm (2X2”), 5X7.5cm (2X3”), 10X10cm (4X4”). **Indications:** WoundClot™ bio-absorbable Surgical Gauze is intended for bleeding control during surgical procedures. **Application:** WoundClot™ Surgical is dissolvable, keep dry until application 1. Open sterile pack. 2. Gently place flat (unfolded) WoundClot™ over the bleeding site. 3. Hold the gauze gently on the bleeding site for a few seconds to achieve adherence, make sure the gauze is saturated with blood (do not move the gauze once placed). 4. The gel formed, will adhere to the wound and stop the bleeding. 5. If the bleeding persists, remove the existing gauze and apply a new gauze over the bleeding site. **Removal:** Although WoundClot™ Surgical is fully absorbable, it is advised to remove excess material once hemostasis is achieved. Saturate injured area with water or saline solution, and gently remove any excess material from the wound. **Storage:** Keep sealed packets in an ambient environment and a dry place for maximum shelf life. **Warnings:** WoundClot™ is not intended for wound healing. WoundClot™ is not antimicrobial or intended to prevent infections. WoundClot™ gauze is supplied sterile. While it is used to help achieve hemostasis in, around, or in proximity to foramina in bone, areas of bony confine, spinal cord, optic nerve and chiasm - excess material must always be removed after hemostasis is achieved since the product could swell and may exert unwanted pressure in these areas.

Do not re-sterilize

WoundClot™ gauze is a single use device. Do not reuse.

Consult instruction for Use

Sterilization using Irradiation

Single sterile barrier system with protective packaging inside

Do not use if package is damaged. Sterility of contents guaranteed unless individual package is damaged or opened.

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CORE SCIENTIFIC CREATIONS

CE 0483

Manufacturer: Core Scientific Creations Ltd.
Office: 29 Nahshon Street, Petah Tikva 4926918, Israel
Manufacturing Site: 6 Bahat Street, Petah Tikva, 4927149, Israel.
T: +972-3-7359555 | F: +972-3-7359544
info@wound-clot.com | www.woundclot.org

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MedNet Swiss GmbH
MedNet EC-REP III GmbH
Produttore: Core Scientific Creations Ltd.