Picture 1

Picture 2

Picture 3

Picture 4

Xmetix LTD.
39 Derech Haifa,
Kiryat Atta, 2822641, Israel

REF TAK7-10-000

60°C (140°F)

30°C (-22°F)

95%

5%

MD

WEBSITE

VIDEO LIBRARY

EN TAK 710 - Instructions for Use

One- or Two-Handed Application
Recommended to place the tourniquet at the topmost part of the limb* (picture 1)
Attache the buckle and firmly pull the strap to tighten it. (picture 2)
Strap should be tight enough that tips of three (3) fingers cannot be slid between the strap and the limb. If the tips of three (3) fingers slide under strap - **retighten**.
Pull the red safety pin out. (picture 3)
Screen displays -"Press ⏻"
Press the Power button. (picture 4)
Screen displays - dynamic bar chart of the tourniquet tightening process, replaced by time counter when tightening phase ends.

Safety Information

Caution
1. To deactivate/activate the backlight, Press the power button.
2. To release or increase pressure on the limb, use the + / - buttons located on the side of the device.

Notes
1. Once applied, TAK 710 should be removed ONLY by medical personnel
2. Dispose the device according to local biohazard procedure.

Limitation of Liability
Xmetix Ltd., its employees, agents, contractors, suppliers, and distributors, assume no liability for injury or damages arising from the application and use of the TAK710. The user assumes all risks of liability.
The TAK710 should only be used as directed by user's military service component guidelines, EMS authority, or under the supervision of a physician.

Standard Use Guidelines

When applied in accordance with directions, the TAK710 is a safe and effective device for controlling life-threatening extremity bleeding.

*If you cannot be sure or cannot take the additional time to examine where the bleeding is coming from based on the situation, the TAK710 can be effectively applied over clothing as high on the arm or leg as possible. The TAK710 must NOT be applied over solid objects within the clothing. As soon as the situation permits, the injured limb should be evaluated and the TAK710 re-positioned 5cm -7.5cm above the injury directly to the skin.

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